
### **MECHANICAL HOIST SLINGS**

## AMPUTEE SLING 22021–22121

#### **MEASUREMENTS**

<b>Shoulder Straps</b> — Leng from material of sling	th 150	250	350
Leg Straps — Leng from material of sling	gth 250	400	550



	Small	Medium	Large	Extra Large
Α	75	80	90	95
В	60	65	70	75
С	70	75	80	85
D	55	65	75	85
E	26	26	26	26
F	125	130	140	145

# AMPUTEE SLING WITH HEAD SUPPORT 22521-22221

Shoulder Straps— Length from material of sling	150	250	350		
Leg Straps — Length from material of sling	250	400	550		
	Small	Medium		Large	Extra Large
Α	85	90	)	105	115
В	60	60 65		70	75
С	70	70 75		80	85
D	55	65		75	85
E	26	20	5	26	26
F	135	14	0	155	165



**IMPORTANT NOTICE TO PURCHASER:** The following is made in lieu of all warranties, expressed or implied. Seller's and manufacturers only obligation shall be to replace such quantity of the product proved to be defective. Neither seller nor manufacturer shall be liable for any injury, loss or damage, direct, indirect or consequential arising out of the inability to use this product. Before using, user shall determine the suitability of the product for the intended use, and user assumes all risk and liability whatsoever in connection therewith. No statement or recommendation not contained herein shall have any force affect unless in agreement signed by officers of seller and manufacturer. All Tech-Assist Products are designed to assist nurses and patients. The patients mursing care should never be neglected, and when Tech-Assist Products are in use, continual assessment should be carried out to ensure the products are still appropriate for the patients condition. Regular maintenance inspection is necessary, and damaged items should be withdrawn from use until rectified.

Unit 2/5 Keele Place KIDMAN PARK, SA 5025 Ph 08 82352200 Fax 08 83533200 sales@tech-assist.com.au

### **TECH-ASSIST PRODUCTS**

www.tech-assist.com.au